

## Obvyklé derivatizační postupy pro jednotlivé skupiny látek

	GC	HPLC, CE, TLC
Aflatoxins	Use HPLC or TLC	Post-column derivatization
Amino acids	n-Butyl ester, TFA	OPA, Dansyl
Aminoglycosides (Antibiotics)	Use HPLC or TLC	FDNB, OPA
Bile acids	Me ester-TMS	p-Bromophenacyl esters
Biogenic amines	Perfluoroacyl	Not usually necessary
Carbohydrates	TFA, TMS, boronate	Post-column
Drugs	Multifunctional.	Phenyl isocyanate
Fatty acids	Me ester	Post-column
Keto acids	Oxime-TMS	Often not necessary
Mycotoxins	See Aflatoxins	
Nitrosamines	Chromatograph without derivatization	Chromatograph without derivatization
Nucleotides, nucleosides and bases	TMS	Depends on specific analyte
Organic acids	TMS	p-Bromophenacyl ester
keto	Oxime-TMS	2,4-DNP
$\alpha$ -keto	Quinoxalinol-silyl	Quinoxalinol
Peptides	Me ester-TFA	Dansyl
Polychlorinated biphenyls (PCB)	No derivatizable groups	No derivatizable groups
Polyamines	TFA	Benzoyl, dansyl
Polycyclic aromatic hydrocarbons (PAH)	(TMS for metabolite-OH groups)	Chromatograph without derivatization
Porphyrins	TMS	Chromatograph without derivatization
Prostaglandins (thromboxanes)	Me ester-oxime-TMS	p-Bromophenacyl ester
Steroids	PFB ester-oxime-TBDMS	
Sugars	MO-TMS	Keto: 2,4-DNP, dansyl
	HFB	benzoyl
	Reduce to alditoles, then	Phenyl isocyanate
	TFA or TMS	Post-column
Vitamins:		
A (retinol)	TMS	Chromatograph without derivatization
B <sub>1</sub> (thiamine)	Acyl	Chromatograph without derivatization, or post-column
B <sub>6</sub>	Use HPLC	Chromatograph without derivatization, or post-column
C	TMS	Chromatograph without derivatization, or oxidize to dehydro-form and prepare 2,4-DNP or quinoxalinol derivatives
D and metabolites	TMS, boronate-TMS	Chromatograph without derivatization
E (tocopherol)	TMS	Chromatograph without derivatization
K	Chromatograph without derivatization (HPLC preferred)	Chromatograph without derivatization

OPA – o-ftalaldehyd

TBDMS – t.-butyldimethylsilyl

FDNB – 1-fluoro-2,4-dinitrobenzen

TFA – trifluoracetyl

TMS – trimethylsilyl

DNP – dinitrofenyl

PFB – pentafluorobenzyl